


% Bo Chen
% 10190141
% 14bc57

% CISC 330
% December 22nd, 2017
% Assignment 4: Gamma Knife Radiosurgery


1. Compute Dose Box



The dose box is the smallest box possible with holds both the PTV and OAR. The functions takes the PTV and the OAR, and returns the smallest dose box defined with the lower left and upper right.



2. Draw 3D Scene

The 3D scene is a visual of the Gamma-Knife radiosurgery. The ellipsoid is the head of the patient, whole the PVT is the red sphere and the OAR is the green sphere. The dose box from question 1 is in red.





3. Compute Linear Function

Taking two 2D points and an x value, the y value is determined. We can use the two points to create a line, and calculate the slope (m) using the point slope formula. This will allow us to solve for the y value.

```
>> Compute_Linear_Function_Test
1.0000

2.0000

3.0000

4.0000

5.0000

6.0000
```

P1 = [1,1] P2 = [5,5] x = {1,2,3,4,6}

The two points $y = x$, so the slope (m) = 1. This means that any x input, should be the output/

4. Compute Radial Distance

Taking a point on a line, the direction vector and a point in space, the distance is found by finding the direction vector from the point on the line and the point in space. We find the shortest distance by taking the cross product of the two direction vectors, and the norm to find the scalar distance.

```
>> Compute_Radial_Distance_Test
Test 1

distance =

    1

Test 2

distance2 =

    1
```

Test 1, The line is vertical, so it should be the difference in x-values


Test 2, The point and line create a right triangle

5. Compute Depth Dose

Using the graph given in the assignment, we can see that at 20mm, the dose is full, while at 120mm the dose is half. Between 0 and d_0 we can use question 3, and similarly for d_0 and beyond.

Depth dose, 20mm = 1, 120mm = 0.5

Depth Dose Function




0	0
1.0000	0.0500
2.0000	1.0900
3.0000	1.0850
4.0000	1.0800
5.0000	1.0750
6.0000	1.0700
7.0000	1.0650
8.0000	1.0600
9.0000	1.0550
10.0000	1.0500
11.0000	1.0450
12.0000	1.0400
.	.
.	.
.	.
110.0000	0.5500
111.0000	0.5450
112.0000	0.5400
113.0000	0.5350
114.0000	0.5300
115.0000	0.5250
116.0000	0.5200
117.0000	0.5150
118.0000	0.5100
119.0000	0.5050
120.0000	0.5000

6. Compute Radial Dose

Similarly to question 5, using the graph given in the question, we can break up the graph into 3 segment and use question 3 at each segment.

Depth dose, $-22.5\text{mm} = 0$, $-7.5 = 1$, $22.5\text{mm} = 0$, $7.5 = 1$



```
>> Compute_Radial_Dose_Test
-23.0000    0
-22.0000    0.0333
-21.0000    0.1000
-20.0000    0.1667
-19.0000    0.2333
-18.0000    0.3000
-17.0000    0.3667
-16.0000    0.4333
-15.0000    0.5000
-14.0000    0.5667
-13.0000    0.6333
-12.0000    0.7000
-11.0000    0.7667
-10.0000    0.8333
-9.0000     0.9000
-8.0000     0.9667
-7.0000     1.0000
-6.0000     1.0000
-5.0000     1.0000
-4.0000     1.0000
-3.0000     1.0000
-2.0000     1.0000
-1.0000     1.0000
0           1.0000
1.0000     1.0000
2.0000     1.0000
3.0000     1.0000
4.0000     1.0000
5.0000     1.0000
6.0000     1.0000
7.0000     1.0000
8.0000     0.9667
9.0000     0.9000
10.0000    0.8333
11.0000    0.7667
12.0000    0.7000
13.0000    0.6333
14.0000    0.5667
15.0000    0.5000
16.0000    0.4333
17.0000    0.3667
18.0000    0.3000
19.0000    0.2333
20.0000    0.1667
21.0000    0.1000
22.0000    0.0333
23.0000    0
```

With the source longitude and latitude on the gamma knife, we can determine the direction vector of the beam. Using azimuth(longitude) and elevation(latitude) we can find points on the sphere. We can use those points subtract the iso center $[0,0,0]$ to find the direction vector.

Using the information from question 7, we can solve for p_1 = iso center - direction vector and find the intersection of the head and the line from the two points.

If the pencil beam intersected or came in contact with the OAR, then the value in the table would become 0 (unsafe).

[illegible]

12. Compute Point Dose from Beam

These number were verified in the Depth & Radial dose tables

```
Computed radiation doses [0,0,0]  
Isocenter point [0,0,0]  
Longitude 0 , Latitude 0
```

```
doseTest1 =
```

```
0.7000
```

```
Longitude 0 , Latitude 90
```

```
doseTest2 =
```

```
0.7000
```

```
Longitude 90 , Latitude 0
```

```
doseTest3 =
```

```
0.6000
```


```
Longitude 90 , Latitude 90
```

```
doseTest4 =
```

```
0.7000
```


14. Compute Dose Volume Histogram

Using doseBox we can take every 2 points (voxel size of 2) and find if they are inside the OAR or PTV. We can use question 12 on all these points to find the radiation dose for the histogram().



15/16. Compute Surface Dose/Compute Dose Surface Histogram

Creating points on a sphere, and calculating the maximum dose in the PTV and OAR.

```
>> Compute_Dose_Surface_Histogram
hottest on PTV
219.6710

33.7500 -2.7245 29.2658

coldest on PTV
156.6270

20.0180 11.0861 16.5679

hottest on OAR
3.3817

12.0688 21.2315 46.5679

coldest on OAR
0

-12.9904 30.0000 37.5000
```

